



630A Great South Road, Greenlane 1051,
P.O. Box 62613, Greenlane 1546
Auckland, New Zealand
Ph: +64 9 580 1959
Fax: +64 9 580 2044,
Email: sales@alphatech.co.nz



Reagent Kit Sets for Wine Analysis



Enzymatic and Colourimetric Kits

Description	Number of tests manual method macro/semi-micro cuvettes	Description	Number of tests manual method macro/semi-micro cuvettes
Acetaldehyde	45 / 135	Anthocyanins	33/110
Acetic Acid	45 / 135	Copper	- / 50
L-ascorbic acid	-/50	Iron	60 / 120
Ammoniacal Nitrogen	50 / 125	Primary Amino Nitrogen	37 /111
Citric Acid	45/135	Total SO ₂	250 / 750
D-gluconic acid	45 / 135	Free SO ₂	40 / 120
D-Glucose/D-Fructose	45 / 135	Tartaric Acid	44 / 132
Glycerol	45 / 135		
L-Lactic Acid	45 / 135		
L-Malic Acid	45 / 135		
Sucrose	-/60		

All kits include comprehensive method information and instructions.

Further information and analyte kits may be found on the following website: www.oenolab.com

See below for spectrophotometer products

See below for Chemetrics Sulfite reagent kit set – no Spectro required.

Please contact Alphatech Systems for all you Wine testing requirements:

Customer Support: sales@alphatech.co.nz Phone (09) 580 1959 or

Graeme Schofield (graeme.schofield@alphatech.co.nz) or Peter Hassan (ph@alphatech.co.nz)



630A Great South Road, Greenlane 1051,
P.O. Box 62613, Greenlane 1546
Auckland, New Zealand
Ph: +64 9 580 1959
Fax: +64 9 580 2044,
Email: sales@alphatech.co.nz

We also offer



Qualitative Reagents for Sulfite wine analysis



Range: 10-100 ppm as SO₂
MDL: 10 ppm / Method: Ripper

Sulfite in Wine Titrets Kit

Cat#

K-9610W

Increments:

10, 11, 12, 13, 14, 15, 16, 18, 20, 25, 30, 35, 40, 50, 70, 100 ppm

Kit comes in a cardboard box and contains everything needed to perform 10 tests: ten ampoules, ten valve assemblies, instructions, and MSDS.



Spectrophotometer Libra S12 UV/VIS Suitable for Wineries

- Single Beam UV/Vis spectrophotometer covering 200-999nm (5nm bandwidth)
- Lamp saver mode
- Customizable user interface and method storage
- Automatic system calibration
- Wide range of sample handling accessories



Further information may be found on the following website : www.biochrom.co.uk

Please contact Alphatech Systems for all you Wine testing requirements:

Customer Support: sales@alphatech.co.nz Phone (09) 580 1959 or

Graeme Schofield (graeme.schofield@alphatech.co.nz) or Peter Hassan (ph@alphatech.co.nz)