

GYROZEN

FASTCOOL EVENT

Gyrozen Special Offer Valid Through Sep. 30, 2013

Save Your Time with
FAST COOL Performance



Buy **1730R** ❄️ Get 24-Hole Fixed Angle Rotor Free!



Refrigerated High Speed
Microcentrifuge

FAST COOL Event Package

FREE



24-Hole Fixed Angle Rotor

Autoclavable, Scratch Resistance,
Chemical Resistance & O-ring Sealed Lid



Tube Adaptor, 0.2 / 0.5 ml

Ordering Information

Cat. No.	Description
GZ-1730R-SPRF-24	1730R, Microcentrifuge with A Fixed Angle Microtube Rotor & Aluminum Lid (24 x 1.5/2 ml) – including 24 ea of 0.2/0.5 ml tube adaptors

GYROZEN

Alphatech Systems Ltd
Phone: 09 580 1959
Email: sales@alphatech.co.nz

www.gyrozen.com

MicroCentrifuge

1730R ❄️



FAST Cool Function



Accurate Speed Control

Max. 17,000 RPM with 24 hole micro rotor



9 Acc / 10 Dec ramps



Comfort Soft Touch Button



Quiet Performance



Up to 100 Programs Memory




**Automatically Recognize
Max. RPM of Loading Rotor**




Max. 17,000 RPM / 27,237 xg

Fast Cool Performance



- High capacity, strong compressor ensures fast cooling of chamber and sample.
- Reach 4°C within 5 minutes for fast start up and fast cool down of sample.

Applications



For All Kinds of Molecular Biology Experiments

Ethanol precipitation, Phenol extraction, NA preparation, and Cell-down, etc.

**RNA/DNA or Protein Filtration using Filter Tubes
(Recommended rotor, GRF-M-m2.0-18)**

Rotors for 1730R



GRF-M-m2.0-24

24 x 1.5/2.0 ml
Max. 17,000 rpm/27,237 xg



NEW



GRF-M-m2.0-18

18 x 2.0 ml/microfilter tube
Max. 15,000 rpm/18,891 xg



GRF-M-m2.0-30

24 x 1.5/2.0 ml
Max. 14,000 rpm/21,124 xg



GRF-M-s0.2-64

64 x 0.2 ml
Max. 12,500 rpm/11,722 xg



GYROZEN www.gyrozen.com

Alphatech Systems Ltd
Phone: 09 580 1959
Email: sales@alphatech.co.nz